



DRYING APPARATUS




08-1000: DRYING APPARATUS, ABDERRHALDEN - For use in vacuum/vapor type drying operations. Consists of one West type condenser, one jacketed drying chamber with S 34/45 joint, one desiccant tube with S 34/45 joint and vacuum connection with one solvent flask.

Code #	Description
01	Complete Assembly
02	Condenser with S 24/40 bottom joint
03	Drying Chamber
04	Desiccant Tube
05	250ml Flask with S 24/40



08-1020: DRYING TUBE, U-SHAPED - Rubber stopper has bent outlet tube. Ends are standard taper inner joints and beaded rim with rubber stopper insert.

Code #	S Joint	S Joint
01	24/40	
02		29/42



08-1040: DRYING TUBE, PLAIN - With bottom hose connection and top tooled for stopper.


Code #	Body Length mm
01	100
02	150
03	200

To Order, Call Specialty Glass, Inc. Toll Free 1-800-899-2137

DRYING APPARATUS

08-1060: DESSICANT COLUMN - Flat o-ring joints on either end of column allow for easy installation and removal of desiccant. Supplied with upper and lower right angle adapters and two size 25 o'rings.

Approximate Size: 55mm OD; 85cm Height



To Order, Call Specialty Glass, Inc. Toll Free 1-800-899-2137



We welcome your requests for custom glass. Please contact us via email, fax or phone for glassware quotes.

sales@specialtyglassinc.com

www.specialtyglassinc.com